



## Foot switch for ophthalmology MFS PHACO-MED Material no.: on request


### Product features

- Pedal with proportional output signal (both vertical and horizontal)
- 4 rocker switches with 2 switching functions each
- Programmable braking points for tactile feedback
- Intuitive and highly precise operability
- Wireless signal transmission via steute wireless SW2.4LE-MED (optional)

### Notes

- Photo shows optional accessories

### Dimensions



### Technical data

#### Applied standards

IEC 60601-1, ANSI/AAMI ES60601-1, IEC 60529

#### Enclosure

thermoplastic, shockproof

#### Rocker switch

robust, operating theatre approved silicone

#### Degree of protection

IPX8 (IEC/EN 60529)

#### Switching system

rocker switch: 8 switching functions  
pedal: with proportional output signal, horizontal:  $-10^\circ \dots 0 \dots 10^\circ$  +/-  $2^\circ$ , vertical:  $0 \dots 15^\circ$  +/-  $1^\circ$

#### Mechanical life

Rocker switch: > 250 000 operations

Pedal: > 1 million operations

#### Connection

2 m cable, other lengths optionally available

#### Switching voltage

max. 25 VAC / 60 VDC

#### Switch-on current

max. 1 A

#### Switching capacity

max. 20 VA

Errors and omissions excepted.